

## RESULTS YOU CAN SEE

BEFORE



Courtesy of Roy Geronemus, MD

AFTER



AFTER 2 SESSIONS



Courtesy of Steve Weiner, MD



AFTER 2 SESSIONS



Courtesy of Guy Cappuccino, MD



AFTER 3 SESSIONS

"I have received a ton of compliments...  
My face and jawline improved dramatically,  
and I am extremely happy with my results!"

— PATIENT G.P.

RESULTS MAY VARY.

## { ANSWERS TO YOUR questions }

### How does the treatment feel?

You will feel some pressure and heat. The treatment is typically tolerated with topical anesthetics. Your practitioner will work with you to find and improve your comfort level.

### How long is a treatment and its downtime?

Depending on the area treated, it can take as little as 15 minutes and up to an hour. There is an initial pinkness and swelling which typically subsides after 24-36 hours. It may take up to a week for this to completely resolve.

### How long before I see results?

Results may vary, but you should see some changes soon after the procedure. As your body rebuilds collagen, gradual improvements occur over the next six months.

### How many treatments will I need?

Your practitioner will discuss your goals, which will determine the length and number of treatments. Dramatic results may require multiple treatments. Typically 2-3 treatments give optimal results.

SEE GENIUS FOR YOURSELF!  
GENIUSFOR.ME

**LUTRONIC**  
Intelligent Care

©2018 LUTRONIC. ALL RIGHTS RESERVED.  
LUTRONIC, Lutronic Intelligent Care, Lutronic Genius,  
and associated logos are service marks, trademarks,  
and/or registered trademarks of Lutronic Corp. and/  
or its subsidiaries in the USA and/or other countries.  
PAT\_Genius\_BRF\_0001\_v1 USA

# INTELLIGENT ANTI-AGING



THAT'S **GENIUS**<sup>™</sup>

# THE SMARTER WAY TO REJUVENATE AGING SKIN

Finally, a way to look refreshed and renewed without having to worry about downtime!




Genius can treat anywhere on the body. Contact your practitioner to find the right treatment for you.

## Long-lasting results.

Unlike some other procedures, Genius helps stimulate your own body into action. Therefore, results may last for years.

## HOW GENIUS works

Lutronic Genius uses unique and intelligent technology to help protect the top layers of the skin while providing real-time feedback during the procedure. This allows for treatments on all skin types and colors, including tanned skin. By delivering customized energy precisely to the dermis, it creates a micro-injury to jump-start your body's natural healing process for firm and renewed skin.



"Genius allows me to provide patients with precise and customized treatments with less pain. My patients love their dramatic, long-lasting results and the quick recovery makes the procedure easily tolerated. I love that I can treat skin of all colors and target skin anywhere on the body."

— S. WEINER, MD

## WHAT IS GENIUS?

Genius treatments help provide lasting unparalleled outcomes with Intelligent RF (Radio Frequency) technology, delivering thermal energy through micro-needles into the dermis. The goal is to stimulate collagen production for impressive and long-lasting results. Achieve youthful skin without the surgical downtime.

The system senses the quality of your skin, and treatment can be adjusted for a great outcome. Genius is the smarter way to treat your skin, and it's only available from Lutronic.

